

Smart-Flow Fluidic Patent Technology

The creative SMART-FLOW fluidic technology is a simple and efficient system, which makes FinoCount-5 with good reliability and free of maintenance.



Output Parameters-Fino-Count-5


25 Basic parameters			
Abbr.	Parameter Name	Abbr.	Parameter Name
WBC	Mean Corpuscular Volume	MCV	White blood cell count
Neu%	Mean Corpuscular Hemoglobin	MCH	Neutrophils percentage
Lym%	Mean Corpuscular Hemoglobin Conce	MCHC	Lymphocytes percentage
Mon%	Red blood cell distribution width coeff	RDW-CV	Monocytes percentage
Eos%	Red blood cell distribution width stand	RDW-SD	Eosinophils percentage
Bas%	Hematocrit	HCT	Basophils percentage
Neu#	Platelet count	PLT	Neutrophils number
Lym#	Mean platelet volume	MPV	Lymphocytes number
Mon#	Platelet distribution width	PDW	Monocytes number
Eos#	Plateletcrit	PCT	Eosinophils number
Bas#	Large platelet count	P-LCC	Basophils number
RBC	Large platelet percentage	P-LCR	Red blood cell
HGB	Hemoglobin Concentration		

4 Research parameters

Abbr.	Parameter Name
ALY%(RUO)	No of Large Immature Cells
LIC%(RUO)	Large Immature Cells
ALY%(RUO)	Abnormal Lymphocytes Numbers
LIC%(RUO)	Large Immature Cells

4 Charts

Abbr.	Parameter Name
WBC Histogram	White blood cell histogram
RBC Histogram	Red Blood cell histogram
PLT Histogram	Platelet Histogram
Diff Scattergram	Differential Scattergram



Minimum Maintenance

- Only require for weekly and annually based maintenance, which is also an extremely simple and cost effective maintenance list.
- Might be the lowest operation cost as well as maintenance difficulty for distributors.

Ingenuous Internal Structure

- Life time tubing - Do not have to be touched ever
- Modular design - Lower maintenance difficulty
- There is no high- Voltage anywhere in the system
- Liquids are separated from electronics
- Valves are easy to reach

Technical Specification

Principles



Parameters

25 Reportable parameters:
WBC, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV, PLT, MPV, PCT, PDW, P-LCR, P-LCC, NEU%, LYM%, MON%, EOS%, BAS%, NEU#, LYM#, MON#, EOS#, BAS#

1 3D Scattergram
3 Histograms(WBC/BASO, RBC, PLT)

4 Research parameter:
ALY%, ALY#, LIC%, LIC#

Test Mode

- CBC mode, CBC+DIFF mode
- Venous whole blood, Capillary whole blood and Prediluted

Throughput

60 tests/hour

Performance

Parameter	Linearity Range	Carry Over	CV
WBC	0-300x10 ⁹ /L	≤0.5%	≤2.0%
RBC	0-8x10 ¹² /L	≤0.5%	≤1.5%
HGB	0-250g/L	≤0.5%	≤1.5%
PLT	0-3000 x10 ⁹ /L	≤1.0%	≤4.0%

Sample Volume

CBC+DIFF mode : ≤20ul
CBC mode : ≤10ul

Data Memory

Up to 100,000 results(including histogram scattergram, patient information)

Display

14 inch touch screen resolution 1366*768

Interface

1 LAN port, 4 USB ports

Communication

Bi-direction LIS, support HL7 protocol
Internal RFID reader

Printout

Support various external USB printers, printout formats user definable

Size/Weight

L * W * H = 480*375*517(mm)
Weight: 36kg

Power Requirement

a.c.100-240V,50/60Hz

Working Environment

- Temperature: 10-30°C
- Humidity: 20% - 85%
- Air pressure: 70~106kPa
- Working latitude: ≤3500m



FinoCount-5

Auto 5 - Part Hematology Analyzer



To Know More About Our Complete Product Range



PSR GROUP OF COMPANIES



Excl. Mktd. by
Finova Healthcare Pvt.Ltd
10, Savitri Chamber, 3rd Floor, D Block., Central Mkt., Prashant Vihar,
Delhi-110085 Ph. : 011-45710537
Email.: sales@psrbiochem.in, Website : www.finovahealthcare.com



Powered By Innovation,
Delivered In Budget

FinoCount-5

Ingenuity for 5-part

FinoCount-5, the new line of 5-part hematology analyzer, is developed by innovative technology, where ingenuity meets advancement. With the innovative technology, simple operations, fewer reagents, intelligent fluidic system design and better performance, all in a smaller more affordable package that will fit any clinician budget and space.

Principle


Tri-angle laser scatter + flow Cytometry + impedance method for WBC.

The 5 part differentiation of the white blood cell can be precisely done by collecting the optical signal when WBC pass through the laser beam.

The front small-angle optical signal can reflect the information of the cell size.

The front large-angle optical signal can reflect the information of nucleus' structure and complexity.


The side angle optical signal can reflect the information of granularity complexity.



Flow Cytometry By Laser


3D Scattergram

• **Rotatable:** Conveniently to see accurate WBC 5-part differential information from any angles.




Dual Methods For BASO Measurement

The first innovative analyzer combined the optical method of BASO(BASO-O) and impedance method of BASO(BASO-I) together, it brings more reliable and stable measurement of BASO pathologic samples and minimized the analysis failure



Dual methods for BASO measurement

Combined optical and Impedance methods, improved the accuracy and stability of BASO measurement, especially for pathologic samples, minimized the analysis failure



Compact Size, Powerful Functions

Compact design with reagents on board the valuable bench space of small labs.



Support most of the printers



LAN port, 4 USB ports

User Friendly Software

- One menu for most of the daily operations
- Single click for troubleshooting
- Flags of abnormal results information

Sample & Patient Information Test Results Shortcut Menus Graphic Information



- 14-inch, High Resolution:1366
- Capacitive Touch screen, Support With Glove

Flags info Status Info

Low Reagents Consumption

Low reagents consumption

035ml Lyse/test, 21.5ml Diluent/test
Approximately 50% less reagent consumption
Higher efficiency on diluent

Probe cleanser consumption

1.5ml daily

3 Reagents for FinoCount-5

W-61D Diluent	20L
W-61LD Lyse	500 ml
W-61LH Lyse	500 ml
W-5P Probe Cleanser	50 ml




User-Oriented, Humanized Design

- Classic, impressive appearance design
- Real-time alerts, one click to remove error

Instrument Status

Laser Status



Smart shutdown system design Intelligent Switch

- Auto switch off (power disconnected) after performing cleaning procedures, no need waiting to switch off manually



True optical system

Solid and Advanced Technology

