

<b>FINO-ESR Fast Detector</b>	
<b>Features</b>	
Better compatibillity compared with westergren method	
Safe Reliable and airtight operation during the whole process from blood collect	
Designed with spirit level of patent technology, which ensure the detector horizc	
10 channels, synchronous operation allowed, only 30 minutes needed for readin	
Stainless steel soleplate with silcgel cushion on which will avoid any noise or bre:	
<b>Technical Specifications</b>	
Application range :	Measuring of erythrocyte sedimentation
range	1-140 mm/h
Measuring time :	30 Mintutes
Analytical capacity	maimum 20 tests/h
Loading capacity :	maimum 10 samples at a time
Loading pattern :	random
Analysis result:	in Westergren sedimentation rate (mm/h)

ion to test free of biological contamination

ontal and aclinic.

ig easy and shift

akage when the esr tube being placed onto the soieplate