

Finou-200 (Semi Automatic Urine Analyzer)

URS-10H Test Items: Urobilinogen, Billirubin, Ketone, Blood, Protein, Nitrile, Leukocytes, Glucose, Specific C

URS-11H Test Items: Urobilinogen, Billirubin, Ketone, Blood, Protein, Nitrile, Leukocytes, Glucose, Specific C

URS-12H Test Items: Urobilinogen, Billirubin, Ketone, Blood, Protein, Nitrile, Leukocytes, Glucose, Specific C

Specification	
Test Principle:	Photoelectric colorimetric method
Optional :	10, 11 or 12 parameters
Test speed :	60-120 strips/ hour
Wavelength :	470nm, 540mm, 620nm
Memory Size:	1000 results, can query the test result according to ID
Display :	192x64 LCD Screen
Port:	Can be connected to urine sediment analysis system and hospital management system through RS232 port
Record mode:	Thermal Prinitng
Environment condition:	15-30 C: RH<70%
Power suppl;y :	AC220V/110V ± 15%, 50/60Hz
Power:	20VA
Weight	3kg

Gravity, PH.

Gravity, PH, Ascorbic Acid.

Gravity, PH, Ascorbic Acid, Microalbumin.