

Finohit Pipette

Features

The FINOHIT single channel manual pipette is among the most ergonomic manual pipettes available.

PIPETTER FEATURES

- Low Spring forces
- Autoclavable liquid end
- Total hand comfort
- Durable Smooth Operation
- Ergonomic
- A Better Tip Ejector
- Accuracy and precision

Technical Specification

1.1 Fixed Volume Pipette

Cat No.	Volume (ul)	Inacc. %	Impr. %
A3012	50	± 2.4	<11.01
A3013	100	± 1.5	<10.81

1.2 Adjustable Volume Pipette

Cat No.	Volume Range (m)	Increment (ul)	Test Volume (ul)	Inacc (%)	impr (%)
A3034	5-50 ul	0.5 ul	50 ul	± 1.5	<11.51
			25 ul	± 2.8	<12.01
			5 ul	± 6.0	<13.01
A3025	10-100 ul	1 ul	100 ul	± 1.5	<11.01
			50 ul	± 2.2	<11.51
			10 ul	± 3.0	<12.01
A3028	100-1000 ul	5 ul	1000 ul	± 1.2	<10.51
			500 ul	± 1.5	<10.51
			100 ul	± 1.5	<11.01



Tip ID
200 ul
200 ul
1000 ul