

# FINGERTIP PULSE OXIMETER

Model No. APO-30C



FDA

CE 0123



- **Artifact removal and anti-motion**
- **Splash proof and drop resistant**
- **Accurate and sensitive**

For The Ones  
You Care Most...



# FINGERTIP PULSE OXIMETER

Model No. APO-30C



## Features:

- OLED display of SpO2, PR, PI Pulse bar and waveform
- Special design with splash-proof and drop-resistant
- Spot check and continuous measuring mode
- Pulse rate analysis for spot check measurement
- UP to 12 groups SpO2 data storage
- Audible and visual alarm with low battery indication
- Automatic Power On / Off
- 4 Direction Display



PSR GROUP OF COMPANIES

  
HEALTHCARE PVT. LTD.

Excl. Mktd. by :

**FINOVA HEALTHCARE PVT. LTD.**

10, Savitri Chamber, 3rd Floor, D Block., Central Mkt., Prashant Vihar  
Delhi-110085 Ph. : 011-41066800, Email.: [finovahealthcare@gmail.com](mailto:finovahealthcare@gmail.com)  
[www.finovahealthcare.com](http://www.finovahealthcare.com)

