

Features :

Non-contact infraerd temperature measurement.
High temperature sound alarm
Fast temperature measurement
Beautiful and exquisite LCD display 5, easy to carry, simple operation.



Category	Project Name	Explain
Specification type	Product Name	Non-contact infrared thermometer
	Product Number	TM-168B
	Specification	/
	Working Environment	-10-40°C
	Temperature Range	32°C-42.9°C
	Basic Accuracy	±0.3°C
	Test Distance	5CM-10CM
	Response Time	500 Millisecond
	Body Material	ABS
	Data Storage	Store 32 Sets of Measurement



CUSTOMER CARE NO.

9654075003



Telephone No.

011-45710537

ADD.: 10, Savitri Chamber, 2nd Floor, D Block, Central Mkt., Prashant Vihar, Delhi-110085,
Email: customercare@psrbiochem.in, sales@psrbiochem.in